

# GuardianCath™

*Protecting Those Who Cannot Protect Themselves*

Color Coded for Patient SAFETY in Identifying the Catheter  
Designed for the Neonatal Patient

## Umbilical Catheterization Tray

- Comprehensive enough to be used as a minor surgical kit
- Umbilical Catheterization, Chest Tube Insertion, Cutdown Tray and compact enough for Transport Teams

## Catheterization Tray

Part Number	Description	Quantity
NCT-100	Catheterization Tray	4/box

### Kit Contains:

1-CSR Wrap  
2-Trays  
1-Tape Measure  
4-Towel-Cloth  
1-Scissors-Iris  
1-Forcep, Adson 1x2  
1-Forcep, Iris 4"  
1-Forcep, Iris 2 1/2"  
1-Hemostat, Straight  
2-Hemostat, Curved  
1-Needle Holder  
1-Probe, 5 1/2 "  
1-Scalpel, #11  
1-Syringe, 1cc  
1-Syringe, 3cc  
1-Syringe, 6cc  
1-Ampule Saline  
1-Needle, 20 GA x 1-1/2  
1-Needle, Filter  
1-Swabstick, Iodophor-3 pack  
1-Umbilical Tie 12"  
4-Gauze, 2"x 2"  
4-Gauze, 4"x 4"  
2-Sharps Containers-Point-Lok™

Manufactured for Marian Medical, Inc. Louisville, KY 40207

Phone: 1.502.895.6991 Fax: 1.502.721.6128

[www.marianmedicalonline.com](http://www.marianmedicalonline.com)