

SafeFEED™

Protecting Those Who Cannot Protect Themselves
Color Coded for Patient **SAFETY** in Identifying the Catheter
Designed for the Neonatal Patient

SafeFEED™ Enteral Only Extension Sets with **Oral** Adapter

- Designed to minimize the chance of inadvertently connecting an enteral feed through an intravenous route
- **"For Enteral Use Only"** is labeled on a tag for additional safety recognition
- Orange stripe to identify the line as an enteral feed

Enteral Only Extension Sets with Oral Adapter

Part Number	Description	Quantity
SFO-36	36" Extension Set	50/Box
SFO-60	60" Extension Set	50/Box
SFOS-60	60" Extension Set w/Stopcock	50/Box
SFOY-60	60" Extension Set w/Y-Site	50/Box

- ◇ **Sterile - Individually Packaged**
- ◇ **Minibore PVC Tubing (.047)**
- ◇ **DEHP and Latex Free**

Manufactured for Marian Medical, Inc. Louisville, KY 40207
Phone: 1.502.895.6991 Fax: 1.502.721.6128
www.marianmedicalonline.com