

GuardianCatheter™

Protecting Those Who Cannot Protect Themselves

Color Coded for Patient **SAFETY** in Identifying the Catheter
Designed for the Neonatal Patient

Infant Silicone Feeding Tube

- Designed as an open-ended catheter with a radio-opaque stripe to monitor the catheters position
- Rounded distal tip for patient comfort and three side holes for added drainage
- Orange funnel and captive closure plug for **SAFETY** in identifying the catheter as an enteral tube

Infant Silicone Feeding Tube with Oral Dose Hub

| Part Number | Description | Quantity |
|--------------------|--------------------|-----------------|
| FTSEN-3540 | 3.5 Fr 40 cm | 10/Box |
| FTSEN-540 | 5.0 Fr 40 cm | 10/Box |
| FTSEN-6540 | 6.5 Fr 40 cm | 10/Box |
| FTSEN-840 | 8.0 Fr 40 cm | 10/Box |
| FTSEN-560 | 5.0 Fr 60 cm | 10/Box |
| FTSEN-6560 | 6.5 Fr 60 cm | 10/Box |
| FTSEN-860 | 8.0 Fr 60 cm | 10/Box |
| FTSEN-590 | 5.0 Fr 90 cm | 10/Box |
| FTSEN-6590 | 6.5 Fr 90 cm | 10/Box |
| FTSEN-890 | 8.0 Fr 90 cm | 10/Box |

Manufactured for Marian Medical, Inc. Louisville, KY 40207

Phone: 1.502.895.6991 Fax: 1.502.721.6128

www.marianmedicalonline.com