

GuardianPICC™

Protecting Those Who Cannot Protect Themselves
Color Coded for Patient **SAFETY** in Identifying the Catheter
Designed for the Neonatal Patient

GuardianPICC™ Peripherally Inserted Central Catheter

- Designed for the neonatal patient
- The PICC line is 28 cm in length with 25 cm of ink markings
- Short extension tail to help decrease catheter breakage and dressing changes

Product Number	Description	Quantity
PICC-2	2Fr Silicone PICC Catheter	5/Box

GuardianPICC™ PICC Procedure Trays

- Neonatal component friendly
- Silicone tourniquet and measuring tape on top for patient assessment
- Neonatal Tegaderm® and 3M steri-strips
- Cloth surgical towels are packed in the procedure tray

GuardianPICC™ PICC Dressing Change Tray

- Neonatal component friendly

Product Number	Description	Quantity
DCT-100	Neonatal PICC Dressing Change Tray	20/Box

Manufactured for Marian Medical, Inc. Louisville, KY 40207

Phone: 1.502.895.6991 Fax: 1.502.721.6128

www.marianmedicalonline.com