

GuardianCatheter™

Protecting Those Who Cannot Protect Themselves

Color Coded for Patient **SAFETY** in Identifying the Catheter
Designed for the Neonatal Patient

Infant Urinary Drainage Catheters and Collection Kits

- Designed as an open-ended catheter with a radio-opaque stripe to monitor the catheters position
- Rounded distal tip for patient comfort and three side holes for added drainage
- Yellow funnel and captive closure plug for **SAFETY** in identifying the catheter as a Urinary Catheter

Urinary Drainage Catheters

Part Number	Description	Quantity
UDC-3540	3.5 Fr 40 cm Silicone Catheter	10/Box
UDC-540	5.0 Fr 40 cm Silicone Catheter	10/Box
UDC-6540	6.5 Fr 40 cm Silicone Catheter	10/Box
UDC-840	8.0 Fr 40 cm Silicone Catheter	10/Box

Urinary Collection Kits

Part Number	Description	Quantity
UCK	Urinary Collection Kit	5/Box

Kit Contains:

- 1-Fenestrated drape
- 1-Towel
- 1-Burette, closed, 150 ml
- 1-Povidone-Iodine swabsticks
- 1-Surgical Lubricant
- 5-Gauze sponges

Manufactured for Marian Medical, Inc. Louisville, KY 40207

Phone: 1.502.895.6991 Fax: 1.502.721.6128

www.marianmedicalonline.com