



Blood Administration Set

Neonatal/Pediatric

Marian Medical, Inc.

Protecting Those That Cannot Protect Themselves

3600 Chamberlain Lane, Suite 122

Louisville, KY 40241

Phone: 502-425-6363

Fax: 502-425-6543

www.marianmedicalonline.com

Blood Administration Set

Neonatal/Pediatric

Marian Medical's Blood Administration Set is a syringe set designed for the preparation and delivery of small aliquots of whole blood, red cells and platelets with a 150-micron filter. It is a cost effective, needleless system that maximizes the usage of adult size units of blood with sterile connection and single donor exposure for the neonatal/pediatric patient. The system comes with a BD™ syringe that is accurately identified and compatible with all IV pumps.

Item Part Number	Item Description
BASN-30/150	Blood Administration Set-w/30 cc syringe and 150 Micron Filter 25 per Box
BASN-60/150	Blood Administration Set-w/60 cc syringe and 150 Micron Filter 25 per Box
BASNWS-60/150	Blood Administration Set-w/out SPIKE-60 cc syringe and 150 Micron Filter 25 per Box
BASN-NS/150	Blood Administration Set-NO Syringe, 150 Micron Filter 50 per Box



- Non-DEHP, Latex free and BPA free blood delivery system
- Sterile, pyrogen-free blood path
- Compatible with sterile connecting technology
- Provides a closed blood delivery system
- Limits the multiple donor exposure to the neonate/pediatric patient
- 150 micron filter in accordance with AABB blood administration requirements
- Needleless system with protected spike for safe handling

