



**Marian Medical, Inc.**

*Protecting Those Who Cannot Protect Themselves*

## Neonatal/Pediatric Blood Administration Products

**Marian Medical, Inc.**

*Protecting Those Who Cannot Protect Themselves*

3600 Chamberlain Lane, Suite 106

Louisville, KY 40241

Phone: 502-425-6363

Fax: 502-425-6543

[www.marianmedicalonline.com](http://www.marianmedicalonline.com)

# Blood Administration Sets

Our syringe sets are designed for the preparation and delivery of all blood and blood products that require the use of a 150-micron filter.

Our syringe sets are a cost effective, needleless system that maximizes the usage of adult size units of blood with sterile connection and single donor exposure.



Part Number	Administration Sets	Quantity
BASN-30/150	Administration Set-w/Spike, 30cc Syringe, 150 Micron Filter ( <i>BD Syringe</i> )	25/Bx
BASN-50/150	Administration Set-w/Spike, 50cc Syringe, 150 Micron Filter ( <i>BD Syringe</i> )	25/Bx
BASN-60/150	Administration Set-w/Spike, 60cc Syringe, 150 Micron Filter ( <i>Monoject Syringe</i> )	25/Bx
BASN-NS/150	Administration Set-w/Spike, 150 Micron Filter, NO Syringe ( <i>For use with any size/manufactures syringe</i> )	50/Bx
BASNWS-30/150	Administration Set-w/OUT Spike, 30cc Syringe, 150 Micron Filter ( <i>For use w/Sterile Docking, BD Syringe</i> )	25/Bx
BASNWS-50/150	Administration Set-w/OUT Spike, 50cc Syringe, 150 Micron Filter ( <i>For use w/Sterile Docking, BD Syringe</i> )	25/Bx
BASNWS-60/150	Administration Set-w/OUT Spike, 60cc Syringe, 150 Micron Filter ( <i>For use w/Sterile Docking, Monoject Syringe</i> )	25/Bx
BA-1507.13	Administration Set-w/Spike, 150 Micron Filter, Y-Port w/Stop Cock, 30cc Syringe, 60 "Small Bore Tube ( <i>BD Syringe</i> )	25/Bx
BMS-2135.17	Administration Set-w/FLL, 150 Micron Filter, MLL	50/Bx
ETT-NW	Exchange Transfusion Set-NICU	1/EA
MMBC105	Cap-Male/Female Red Luer Lock	100/Bx

*\*BD is a trademark of Becton Dickinson. \*Monoject is a trademark of Covidien/Cardinal Healthcare*