

# ENFit Feeding Tube Accessories

## ENFit Drainage Bags

Enteral drainage bags connect directly to feeding tubes or to ENFit extension sets (connected to feeding devices) to allow collection of drained fluids.



Medicina code	Product Description	Box Qty
MED50L	50ml ENFit Drainage Bag	25
MED100L	100ml ENFit Drainage Bag	25
MED250L	250ml ENFit Drainage Bag	25
MED500L	500ml ENFit Drainage Bag	25
MED1000L	1000ml ENFit Drainage Bag	25
MED51L	50ml ENFit Drainage Bag including 75cm extension tubing & bag strap	20
MED101L	100ml ENFit Drainage Bag including 75cm extension tubing & bag strap	20
MED251L	250ml ENFit Drainage Bag including 75cm extension tubing & bag strap	20
MED501L	500ml ENFit Drainage Bag including 150cm extension tubing & bag strap	20
MED1001L	1000ml ENFit Drainage Bag including 150cm extension tubing & bag strap	20
MED1000W	1000ml Drainage Bag Including 100cm wide bore cone adapter extension tubing & bag strap	20

### Features


- **Sterile**
- **Individually packaged**
- **Purple connector**
- **ENFit connector**
- **Extract samples of gastric fluid for analysis**
- **Clamp to prevent leaking**
- **Antibacterial filter allows gas to escape so the bag collects 100% stomach contents**
- **Three drainage bag options:**
  - **Standard Drainage Bag**
  - **Drainage Bag including extension tubing and bag strap**
  - **Drainage Bag including wide bore cone adapter extension tubing & bag strap (compatible with RT range – page 14)**
- **DEHP/Latex free**



# ENFit Decompression System

The ENFit Gastric Decompression System removes gas from a patient's stomach, eliminating the pain and discomfort associated with fluid and gas retention.

Patients can digest enteral formula at their own rate and reach their calorific goals.



## Features

- Sterile
- Individually packaged
- Purple connector
- ENFit connector
- Relieves gastric distention (bloating of the stomach)
- Eliminates the pain and discomfort associated with fluid and gas retention
- Allows patients to digest enteral formula at their own rate and reach their enteral calorific goals
- Available in two sizes (Children & Adults)
- DEHP/Latex free

Medicina code	Product Description	Box Qty
EV01L	250ml ENFit Gastric Decompression System	30
EV02L	500ml ENFit Gastric Decompression System	30

