

## ENFit Reusable Syringes



### Features

- Reusable/Washable
- Individually packaged
- Purple
- ENFit compliant
- Clear measurement markings
- Low dose tip design to improve dose accuracy (1ml and 2.5ml)
- Three piece syringe (barrel, plunger and silicone gasket)
- O-ring silicone gasket to prevent swelling when cleaning
- DEHP/Latex Free

Medicina code	Product Description	Box Qty
LHE01LD	1ml ENFit Reusable Low dose Syringe	100
LHE25LD	2.5ml ENFit Reusable Low Dose Syringe	100
LHE05	5ml ENFit Reusable Syringe	100
LHE10	10ml ENFit Reusable Syringe	100
LHE20	20ml ENFit Reusable Syringe	80
LHE60	60ml ENFit Reusable Syringe	60

## Oral Tip Low Dose Syringes (ENFit compatible)

The oral tip low dose syringe has been designed to give feed and medication orally, as well as being compatible with ENFit devices.



### Features

- Sterile
- Individually packaged
- Purple
- ENFit compatible
- Clear measurement markings
- Low dose tip design to improve dose accuracy
- Three piece syringe (barrel, plunger and silicone gasket)
- Smooth oral tip connector
- DEHP/Latex Free

Medicina code	Product Description	Box Qty
OT01LD	1ml Oral Tip Low-Dose Syringe (ENFit Compatible)	100
OT25LD	2.5ml Oral Tip Low Dose Syringe (ENFit Compatible)	100