


ENFit Extension Sets

Several long and short ENFit extension sets are available (with and without Y-Sites) to connect an ENFit device/pump to a feeding tube.



Features

- **Sterile**
- **Individually packaged**
- **Purple connectors**
- **ENFit connectors**
- **Extension sets provide longer length tubing to connect to a feeding tube**
- **Wide bore extension sets are available for high density feeds**
- **Y-site option for administration of medication**
- **DEHP/ Latex Free**

Medicina code	Product Description	Box Qty
EP015L	15cm Long Extension Set ENFit	50
EP015LYW	15cm Long Y Site Extension Set ENFit Wide Bore	50
EP075L	75cm Long Extension Set ENFit	50
EP075LY	75cm Long Y Site Extension Set ENFit	50
EP075LW	75cm Long Extension Set ENFit Wide Bore	50
EP075LYW	75cm Long Y Site Extension Set ENFit Wide Bore	50
EP150L	150cm Long Extension Set ENFit	50
EP150LY	150cm Long Y Site Extension Set ENFit	50

