


Nasogastric & Nasojejunal feeding tubes

Medicina has a comprehensive range of nasogastric and nasojejunal tubes to feed a wide variety of patients from neonates to adults both in hospital and at home.

Fully NPSA compliant


In March 2011 the NPSA UK issued a Patient Safety Alert which updated the NPSA Alert from March 2005. These alerts are very clear in the physical characteristics needed for a tube that is going to be used for feeding, flushing or for delivering medication. All Medicina nasogastric and nasojejunal tubes are fully compliant with both alerts.

ENFit Polyurethane Nasogastric Feeding Tubes



Features

- Sterile
- ENFit compliant
- Polyurethane
- Radio-opaque stripe
- Centimetre markings throughout the length of the tube
- Two side eye holes
- Colour coded strap to signify french size of tube
- Tethered cap
- Rounded tip to aid insertion
- Weighted tubes available
- Up to 90 days use

 Hospital guidelines must be followed for nasal erosion and ventilated patients.


- DEHP/Latex Free

Medicina code	Product Description	Box Qty
LG4/40	4FR x 40cm	50
LG4/50	4FR x 50cm	50
LG4/80	4FR x 80cm	50
LG5/40	5FR x 40cm	50
LG5/50	5FR x 50cm	50
LG5/80	5FR x 80cm	50
LG6/40	6FR x 40cm	50
LG6/50	6FR x 50cm	50
LG6/80	6FR x 80cm	50
LG6/90	6FR x 90cm	50
LG8/40	8FR x 40cm	50
LG8/50	8FR x 50cm	50
LG8/80	8FR x 80cm	50
LG8/90	8FR x 90cm	50
LG8/120	8FR x 120cm	50
LG10/80	10FR x 80cm	50
LG10/120	10FR x 120cm	50
LG12/80	12FR x 80cm	50
LG12/120	12FR x 120cm	50
LG14/80	14FR x 80cm	50
LG14/120	14FR x 120cm	50
LG6/50W	6FR x 50cm (weighted)	50
LG6/80W	6FR x 80cm (weighted)	50
LG8/50W	8FR x 50cm (weighted)	50
LG8/80W	8FR x 80cm (weighted)	50
LG10/80W	10FR x 80cm (weighted)	50
LG12/80W	12FR x 80cm (weighted)	50

Additional size tubes are available upon request.



ENFit Polyurethane Nasogastric Feeding Tubes with Guidewire



Features

- Sterile
- ENFit compliant
- Polyurethane
- Fully radiopaque
- Barium sulphate tip – highlights end of tube clearly under X-Ray
- Centimetre markings throughout the length of the tube
- Two side eye holes
- Colour coded strap to identify French size of tube
- Tethered cap
- Rounded tip to aid insertion
- Guidewire to aid insertion (ENFit connection)
- Includes documentation label
- Weighted tubes available
- Up to 90 days use



Hospital guidelines must be followed for nasal erosion and ventilated patients.


- DEHP/Latex Free

UpTo
90

Medicina code	Product Description	Box Qty
NGP6/55L	6FR x 55cm	25
NGP6/75L	6FR x 75cm	25
NGP6/85L	6FR x 85cm	25
NG6/100L	6FR x 100cm	25
NGP8/55L	8FR x 55cm	25
NGP8/75L	8FR x 75cm	25
NGP8/85L	8FR x 85cm	25
NG8/100L	8FR x 100cm	25
NGP8/120L	8FR x 120cm	25
NGP10/85L	10FR x 85cm	25
NG10/100L	10FR x 100cm	25
NGP10/120L	10FR x 120cm	25
NGP12/100L	12FR x 100cm	25
NGP12/120L	12FR x 120cm	25
NGP14/120L	14FR x 120cm	25
NGP8/85LW	8FR x 85cm (weighted)	25
NGP10/85LW	10FR x 85cm (weighted)	25



ENFit Polyurethane Y-Port Nasogastric Feeding Tubes with Guidewire



Features

- Sterile
- Y-Port ENFit compliant
- Polyurethane
- Fully radiopaque
- Barium sulphate tip – highlights end of tube clearly under X-Ray
- Centimetre markings throughout the length of the tube
- Two side eye holes
- Colour coded strap to identify French size of tube
- Tethered cap
- Rounded tip to aid insertion
- Guidewire to aid insertion (ENFit connection)
- Includes documentation label
- Up to 90 days use



Hospital guidelines must be followed for nasal erosion and ventilated patients.

- DEHP/Latex Free

Medicina code	Product Description	Box Qty
NGP5/55LY	5FR x 55cm	25
NGP6/55LY	6FR x 55cm	25
NGP6/75LY	6FR x 75cm	25
NGP6/85LY	6FR x 85cm	25
NGP6/100LY	6FR x 100cm	25
NGP8/55LY	8FR x 55cm	25
NGP8/75LY	8FR x 75cm	25
NGP8/85LY	8FR x 85cm	25
NGP8/100LY	8FR x 100cm	25
NGP8/120LY	8FR x 120cm	25
NGP10/85LY	10FR x 85cm	25
NGP10/100LY	10FR x 100cm	25
NGP10/120LY	10FR x 120cm	25
NGP12/85LY	12FR x 85cm	25
NGP12/100LY	12FR x 100cm	25
NGP12/120LY	12FR x 120cm	25
NGP14/85LY	14FR x 85cm	25
NGP14/120LY	14FR x 120cm	25
NGP16/100LY	16FR x 100cm	25

UpTo
90



ENFit Polyurethane Nasojejunal Feeding Tubes with Guidewire

The nasojejunal tubes have an open end feeding port that allows them to be passed over an endoscopically or radiologically placed guidewire.



Features

- Sterile
- ENFit compliant
- Polyurethane
- Fully radiopaque
- Centimetre markings throughout the length of the tube
- One side eye hole
- Colour coded strap to identify French size of tube
- Tethered cap
- Open tip (Excluding weighted tubes)
- Guidewire to aid insertion (ENFit connector)
- Includes documentation label
- Weighted tubes available
- Up to 90 days use



Hospital guidelines must be followed for nasal erosion and ventilated patients.

- DEHP/Latex Free


UpTo
90



ENFit Feeding/Aspiration Tubes

These tubes are primarily designed to aspirate stomach contents in adults and children.

They can also be used for continuous or intermittent nasogastric enteral feeding, administering medication and flushing.



Features

- Sterile
- ENFit compliant
- Polyurethane
- ENFit and funnel connector
- Radio-opaque stripe
- Centimetre markings throughout the length of the tube
- Four large side eye holes
- Colour coded strap to signify french size of tube
- Tethered cap
- Rounded tip to aid insertion
- Includes documentation label
- Weighted tubes available
- Tube including guidewire available
- Up to 90 days use



Hospital guidelines must be followed for nasal erosion and ventilated patients.

- DEHP/Latex Free

Medicina code	Product Description	Box Qty
RT6/80L	6FR x 80cm	25
RT8/80L	8FR x 80cm	25
RT10/100L	10FR x 100cm	25
RT12/100L	12FR x 100cm	25
RT14/100L	14FR x 100cm	25
RT16/100L	16FR x 100cm	25
RT18/100L	18FR x 100cm	25
RT20/100L	20FR x 100cm	25
RT8/120LM	8FR x 120cm (Including Guidewire)	25
RT10/120LM	10FR x 120cm (Including Guidewire)	25
RT12/120LM	12FR x 120cm (Including Guidewire)	25
RT14/120LM	14FR x 120cm (Including Guidewire)	25
RT16/120LM	16FR x 120cm (Including Guidewire)	25
RT18/120LM	18FR x 120cm (Including Guidewire)	25
RT8/120LMW	8FR x 120cm (Including Guidewire & Weights)	25
RT10/120LMW	10FR x 120cm (Including Guidewire & Weights)	25
RT12/120LMW	12FR x 120cm (Including Guidewire & Weights)	25
RT14/120LMW	14FR x 120cm (Including Guidewire & Weights)	25
RT16/120LMW	16FR x 120cm (Including Guidewire & Weights)	25
RT18/120LMW	18FR x 120cm (Including Guidewire & Weights)	25

UpTo
90



Medicina Specialist Tubes



ENFit Double Lumen Feeding/Aspiration Tubes

The ENFit double lumen nasojejunal feeding and gastric aspiration tube is used to feed into the small bowel and simultaneously aspirate the stomach. There are two versions of the tube depending on the method of placement.

The EF2L has an internal guidewire for surgical or endoscopic placement.


The EF2EL has an external guidewire and an open tip for an "Over the Wire" placement technique either radiologically or by using the guidewire down the centre of the endoscope.

The aspiration lumen of these tubes may be used to feed the patient as they recover. Care must be taken changing the labels on the tubes.



Features

- Sterile
- ENFit compliant
- Double lumen tube
- Polyurethane
- Radio-opaque stripe
- Centimetre markings throughout the length of the tube
- Internal or external guidewire depending on method of placement
- Rounded tip or open ended tip depending on method of placement
- 14FR/ 10FR x 130cm tube. 14FR x 90cm to stomach exit hole (double lumen) and 10FR x 40cm to jejunal exit hole (single lumen).
- Up to 90 days use

 Hospital guidelines must be followed for nasal erosion and ventilated patients.

- DEHP/Latex Free


Medicina code	Product Description	Box Qty
EF2L	14FR/ 10FR x 130cm ENFit Double Lumen Feeding / Aspiration Tube with internal guidewire and closed tip	5
EF2EL	14FR/ 10FR x 130cm ENFit Double Lumen Feeding / Aspiration Tube with external guidewire and an open tip	5

UpTo
90

ENFit Double Lumen Salem Tubes

The “SALEM” type tubes are designed for simultaneous gastric aspiration and ventilation.

The double lumen limits depression of the stomach and protects the stomach lining. They can also be used for continuous or intermittent nasogastric enteral feeding, administering medication and flushing.



Features

- Sterile
- ENFit Compliant
- Double lumen for Aspiration/ Ventilation
- Polyurethane
- Full length radiopaque stripe
- Centimetre markings throughout the length of the tube
- Four large side eye holes
- Rounded tip to aid insertion
- Colour coded to signify french size of tube
- Up to 90 days use



Hospital guidelines must be followed for nasal erosion and ventilated patients.

- DEHP/Latex free

Medicina code	Product Description	Box Qty
25010ML	10FR x 60cm ENFit Salem Tube	10
25012ML	12FR x 120cm ENFit Salem Tube	10
25014ML	14FR x 120cm ENFit Salem Tube	10
25016ML	16FR x 120cm ENFit Salem Tube	10

UpTo
90