

ENFit Gravity Sets

Medicina have a range of ENFit Gravity Feed Sets including a Universal Bottle Adapter or ENFit Cross Spike.



Features

- Sterile
- Individually packaged
- Purple connectors
- ENFit connector
- Available with Universal Bottle Adapter or ENFit Cross Spike
- Connects to ENFit feeding tubes
- Drip Chamber
- ENFit Y Site
- Roller Clamp
- The MGS12L and MGS13L include a detachable Funnel Adapter. The detachable Funnel Adapter can become a part of a gastrostomy feeding device making it ENFit compatible
- DEHP/Latex free

Medicina code	Product Description	Box Qty
MGS10L	ENFit Gravity Set with Universal Bottle Adapter	30
MGS11L	ENFit Gravity Set with ENFit Cross Spike	30
MGS12L	ENFit Gravity Set with Universal Bottle Adapter & Detachable Funnel Adapter	30
MGS13L	ENFit Gravity Set with ENFit Cross Spike & Detachable Funnel Adapter	30

